

# SRA1 Antibody



PACO12452

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SRA1

**Uniprot**

Q9HD15

**Synonyms:**

steroid receptor RNA activator 1; SRA1; MGC87674; SRA; SRAP; STRAA1; pp7684 ;

**Immunogen:**

Human SRA1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A