

# SRGAP3 Antibody



PACO12458

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SRGAP3

**Uniprot**

O43295

**Synonyms:**

SLIT-ROBO Rho GTPase activating protein 3; SRGAP3; ARHGAP14; KIAA0411; MEGAP; SRGAP2; WRP ;

**Immunogen:**

Human SRGAP3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A