

SRRM1 Antibody



PACO12469

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SRRM1

Uniprot

Q8IYB3

Synonyms:

serine/arginine repetitive matrix 1; SRRM1;160-KD; MGC39488; POP101; SRM160 ;

Immunogen:

Human SRRM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A