

SSH3 Antibody



PACO12475

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SSH3

Uniprot

Q8TE77

Synonyms:

slingshot homolog 3 (Drosophila); SSH3; FLJ10928; FLJ20515; SSH-3 ;

Immunogen:

Human SSH3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A