

SSR4 Antibody



PACO12478

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SSR4

Uniprot

P51571

Synonyms:

signal sequence receptor, delta (translocon-associated protein delta); SSR4; TRAPD ;

Immunogen:

Human SSR4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A