

SSX5 Antibody



PACO12485

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SSX5

Uniprot

O60225

Synonyms:

synovial sarcoma, X breakpoint 5; SSX5; MGC9494 ;

Immunogen:

Human SSX5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A