

ST7L Antibody



PACO12495

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ST7L

Uniprot

Q8TDW4

Synonyms:

suppression of tumorigenicity 7 like; ST7L; FAM4B; FLJ11657; FLJ20284; ST7R; STLR ;

Immunogen:

Human ST7L

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A