

# ST8SIA3 Antibody



PACO12496

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ST8SIA3

**Uniprot**

O43173

**Synonyms:**

ST8 alpha-N-acetyl-neuraminide alpha-2,8-sialyltransferase 3; ST8SIA3; SIAT8C; ST8SialII ;

**Immunogen:**

Human ST8SIA3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A