

STAG2 Antibody



PACO12500

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STAG2

Uniprot

Q8N3U4

Synonyms:

stromal antigen 2; STAG2; DKFZp686P168; DKFZp781H1753; FLJ25871; SA-2; SA2; bA517O1.1 ;

Immunogen:

Human STAG2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A