

STARD8 Antibody



PACO12509

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STARD8

Uniprot

Q92502

Synonyms:

StAR-related lipid transfer (START) domain containing 8; STARD8; DKFZp686H1668; DLC3; KIAA0189; STARTGAP3 ;

Immunogen:

Human STARD8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A