

STIM2 Antibody



PACO12523

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STIM2

Uniprot

Q9P246

Synonyms:

stromal interaction molecule 2; STIM2; FLJ39527; KIAA1482 ;

Immunogen:

Human STIM2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A