

STK17A Antibody



PACO12529

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STK17A

Uniprot

Q9UEE5

Synonyms:

serine/threonine kinase 17a; STK17A; DRAK1 ;

Immunogen:

Human STK17A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A