

STK32C Antibody



PACO12534

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STK32C

Uniprot

Q86UX6

Synonyms:

serine/threonine kinase 32C; STK32C; MGC23665; PKE; RP11-140A10.1; YANK3 ;

Immunogen:

Human STK32C.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A