

STOML3 Antibody



PACO12545

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STOML3

Uniprot

Q8TAV4

Synonyms:

stomatin (EPB72)-like 3; STOML3; Epb7.2l; SRO ; SLP-3, SLP3

Immunogen:

Human STOML3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A