

STRN Antibody



PACO12549

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STRN

Uniprot

O43815

Synonyms:

striatin, calmodulin binding protein; STRN; MGC125642; SG2NA ;

Immunogen:

Human STRN.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A