

STX10 Antibody



PACO12553

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

STX10

Uniprot

O60499

Synonyms:

syntaxin 10; STX10; SYN10; hsyn10 ;

Immunogen:

Human STX10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A