

STX12 Antibody



PACO12555

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

STX12

Uniprot

Q86Y82

Synonyms:

syntaxin 12; STX12; MGC51957; STX13; STX14 ;

Immunogen:

Human STX12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A