

STYK1 Antibody



PACO12569

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

STYK1

Uniprot

Q6J9G0

Synonyms:

serine/threonine/tyrosine kinase 1; STYK1; DKFZp761P1010; NOK; SuRTK106 ;

Immunogen:

Human STYK1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A