

# SULT1A1 Antibody



PACO12576

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SULT1A1

**Uniprot**

P50225

**Synonyms:**

sulfotransferase family, cytosolic, 1A, phenol-preferring, member 1; SULT1A1; HAST1/HAST2; MGC131921; MGC5163; P-PST; PST; ST1A3; STP; STP1; TSPST1 ;

**Immunogen:**

Human SULT1A1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A