

# SULT1A3 Antibody



PACO12577

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SULT1A3

**Uniprot**

P50224

**Synonyms:**

sulfotransferase family, cytosolic, 1A, phenol-preferring, member 3; SULT1A3; HAST; HAST3; M-PST; MGC117469; ST1A5; STM; SULT1A4; TL-PST ;

**Immunogen:**

Human SULT1A3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A