

SULT1B1 Antibody



PACO12578

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SULT1B1

Uniprot

O43704

Synonyms:

sulfotransferase family, cytosolic, 1B, member 1; SULT1B1; MGC13356; ST1B2; SULT1B2 ;

Immunogen:

Human SULT1B1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A