

SULT2B1 Antibody



PACO12582

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SULT2B1

Uniprot

O00204

Synonyms:

sulfotransferase family, cytosolic, 2B, member 1; SULT2B1; HSST2 ;

Immunogen:

Human SULT2B1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A