

SULT4A1 Antibody



PACO12583

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SULT4A1

Uniprot

Q9BR01

Synonyms:

sulfotransferase family 4A, member 1; SULT4A1; BR-STL-1; BRSTL1; DJ388M5.3; MGC40032; NST; SULTX3; hBR-STL-1 ;

Immunogen:

Human SULT4A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A