

SUMF2 Antibody



PACO12584

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SUMF2

Uniprot

Q8NBJ7

Synonyms:

sulfatase modifying factor 2; SUMF2; DKFZp566l1024; DKFZp686l1024; DKFZp686L17160; DKFZp781L1035; MGC99485; pFGE ;

Immunogen:

Human SUMF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A