

SURF2 Antibody



PACO12589

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SURF2

Uniprot

Q15527

Synonyms:

surfeit 2; SURF2; SURF-2 ;

Immunogen:

Human SURF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A