

# TAF1L Antibody



PACO12626

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TAF1L

**Uniprot**

Q8IZX4

**Synonyms:**

TAF1 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 210kDa-like; TAF1L; MGC134910; TAF2A2 ;

**Immunogen:**

Human TAF1L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A