

TAF7 Antibody



PACO12629

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TAF7

Uniprot

Q15545

Synonyms:

TAF7 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 55kDa; TAF7; TAF2F; TAFII55 ;

Immunogen:

Human TAF7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A