

TAOK2 Antibody



PACO12636

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TAOK2

Uniprot

Q9UL54

Synonyms:

TAO kinase 2; TAOK2; KIAA0881; MAP3K17; PSK; PSK1; TAO1; TAO2 ;

Immunogen:

Human TAOK2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A