

TAOK3 Antibody



PACO12637

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TAOK3

Uniprot

Q9H2K8

Synonyms:

TAO kinase 3; TAOK3; DKFZp666H245; DPK; FLJ31808; JIK; MAP3K18 ;

Immunogen:

Human TAOK3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A