

TASP1 Antibody



PACO12641

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TASP1

Uniprot

Q9H6P5

Synonyms:

taspase, threonine aspartase, 1; TASP1; C20orf13; FLJ20212; MGC39159; dJ585114.2 ;

Immunogen:

Human TASP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A