

TAX1BP1 Antibody



PACO12642

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TAX1BP1

Uniprot

Q86VP1

Synonyms:

Tax1 (human T-cell leukemia virus type I) binding protein 1; TAX1BP1; CALCOCO3; T6BP; TXBP151 ;

Immunogen:

Human TAX1BP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A