

TBC1D22A Antibody



PACO12644

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TBC1D22A

Uniprot

Q8WUA7

Synonyms:

TBC1 domain family, member 22A; TBC1D22A; C22orf4; HSC79E021 ;

Immunogen:

Human TBC1D22A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A