

TBC1D23 Antibody



PACO12645

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TBC1D23

Uniprot

Q9NUY8

Synonyms:

TBC1 domain family, member 23; TBC1D23; DKFZp667G062; FLJ11046; NS4ATP1 ;

Immunogen:

Human TBC1D23.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A