

# TBL1X Antibody



PACO12652

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TBL1X

**Uniprot**

O60907

**Synonyms:**

transducin (beta)-like 1X-linked; TBL1X; EBI; SMAP55; TBL1 ;

**Immunogen:**

Human TBL1X.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A