

TBRG4 Antibody



PACO12658

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

TBRG4

Uniprot

Q969Z0

Synonyms:

transforming growth factor beta regulator 4; TBRG4; CPR2; FASTKD4; KIAA0948 ;

Immunogen:

Human TBRG4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A