

TCEA3 Antibody



PACO12665

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TCEA3

Uniprot

O75764

Synonyms:

transcription elongation factor A (SII), 3; TCEA3; TFIS; TFIS. H ;

Immunogen:

Human TCEA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A