

TCEB2 Antibody



PACO12669

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ELOB

Uniprot

Q15370

Synonyms:

transcription elongation factor B (SIII), polypeptide 2 (18kDa, elongin B); TCEB2; ELOB; SIII ;

Immunogen:

Human TCEB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A