

TCF12 Antibody



PACO12671

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TCF12

Uniprot

Q99081

Synonyms:

transcription factor 12; TCF12; HEB; HTF4; HsT17266; bHLHb20 ;

Immunogen:

Human TCF12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A