

# TCF7L2 Antibody



PACO12674

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TCF7L2

**Uniprot**

Q9NQB0

**Synonyms:**

transcription factor 7-like 2 (T-cell specific, HMG-box); TCF7L2; TCF-4; TCF4 ;

**Immunogen:**

Human TCF7L2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A