

TCOF1 Antibody



PACO12678

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TCOF1

Uniprot

Q13428

Synonyms:

Treacher Collins-Franceschetti syndrome 1; TCOF1; MFD1; treacle ;

Immunogen:

Human TCOF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A