

TCP1 Antibody



PACO12679

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TCP1

Uniprot

P17987

Synonyms:

t-complex 1; TCP1; CCT-alpha; CCT1; CCTa; D6S230E; TCP-1-alpha ;

Immunogen:

Human TCP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A