

TCP10L Antibody



PACO12681

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TCP10L

Uniprot

Q8TDR4

Synonyms:

t-complex 10 (mouse)-like; TCP10L; C21orf77; FLJ10932; MGC34297; PRED77; TCP10A-2 ;

Immunogen:

Human TCP10L.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A