

TCTN3 Antibody



PACO12685

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF, IP

Recommended dilutions:

Protein Background:

Gene ID:

TCTN3

Uniprot

Q6NUS6

Synonyms:

tectonic family member 3; TCTN3; C10orf61; DKFZp564D116; TECT3 ; Tectonic-3;

Immunogen:

Human TCTN3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A