

TEAD1 Antibody



PACO12690

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TEAD1

Uniprot

P28347

Synonyms:

TEA domain family member 1 (SV40 transcriptional enhancer factor); TEAD1; AA; REF1; TCF13; TEF-1 ;

Immunogen:

Human TEAD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A