

TERF1 Antibody



PACO12695

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TERF1

Uniprot

P54274

Synonyms:

telomeric repeat binding factor (NIMA-interacting) 1; TERF1; FLJ41416; PIN2; TRBF1; TRF; TRF1; hTRF1-AS; t-TRF1 ;

Immunogen:

Human TERF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A