

# TEX10 Antibody



PACO12701

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TEX10

**Uniprot**

Q9NXF1

**Synonyms:**

testis expressed 10; TEX10; FLJ20287; bA208F1.2 ;

**Immunogen:**

Human TEX10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A