

# TEX13A Antibody



PACO12704

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

TEX13A

**Uniprot**

Q9BXU3

**Synonyms:**

testis expressed 13A; TEX13A; MGC131984 ;

**Immunogen:**

Human TEX13A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A