

PACO12723

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TGDS

Uniprot

O95455

Synonyms:

TDP-glucose 4,6-dehydratase; TGDS; SDR2E1; TDPGD ;

Immunogen:

Human TGDS.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A