

TGFB1 Antibody



PACO12724

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TGFB1

Uniprot

P01137

Synonyms:

Transforming growth factor beta-1, TGF-beta-1, LAP, TGFB1, TGFB

Immunogen:

Human TGFB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A