

TGFB3 Antibody



PACO12726

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TGFB3

Uniprot

P10600

Synonyms:

transforming growth factor, beta 3; TGFB3; ARVD; FLJ16571; TGF-beta3 ;

Immunogen:

Human TGFB3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A